# System II Hex Cassette

EAD-0102, EAD-0104, EAD-0105, EAD-0106, EAD-0107,  
EAD-0108, EAD-0109, EAD-0110, EAD-0111, EAD-0112,  
EAD-0115, EAD-0116, EAD-0118, EAD-0119, EAD-0120

Issue Date: FEB 2023 01.087.863 V4

EMERGO EUROPE  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands

CellPath Ltd  
80 Mochdre Enterprise Park,  
Newtown, Powys, SY16 4LE, UK  
T: +44 (0)1686 611333  
E: sales@cellpath.com  
cellpath.com

# System II Hex Cassette

**EAD-0102, EAD-0104, EAD-0105, EAD-0106, EAD-0107, EAD-0108, EAD-0109, EAD-0110, EAD-0111, EAD-0112, EAD-0115, EAD-0116, EAD-0118, EAD-0119, EAD-0120**

A small device intended to be used in the laboratory to support / contain clinical tissue specimens (e.g., biopsy samples) to facilitate specimen processing [e.g., fixation (encapsulation in paraffin wax), dehydration, infiltration] in preparation for subsequent cytological or histological examination, and storage. It is in the form of a synthetic polymer container, with or without a lid, with drainage holes to maximize fluid exposure during tissue submersion in a parent device into which the container is loaded. This is a single-use device.

For laboratory use.

Store in a dry place. Don't store in direct sunlight.

This product does not expire under normal environmental conditions.

Marking capabilities: EpreDia Printmate, Leica IPC, Sakura Autowrite, CellMark pens, 2B pencils.

Processing Limitations: Suitable for high temperature processing. Suitable for microwave processing.

## Instructions for Use

1. Insert the hopper tube into the EpreDia PrintMate cassette printer.

1. Insert the thread into the Leica IPC cassette printer. Holding the top cassette pull the thread up to release.

1. Insert cassettes manually into printer loading system.

2. Place tissue into the System II Cassette. We recommend that only one sample grossed between 2mm and 4.5mm thick is placed in the cassette.

3. Using either a plastic or metal lid, attach the lid. Close the lid and place the cassette in a fixative ensuring it is fully submerged.

NOTE: Compatible lids are CellPath plastic System II Hex Cassette Lids [EAE-0102-03A] or CellPath Stainless Steel Lids [EAW-0114-08A].

4. Process the cassette through solvents according to your usual protocol until the specimens are ready for embedding. If you are using any unusual reagents, we suggest that you carry out preliminary chemical compatibility tests before using System II Cassette to contain tissue samples.

4. Carefully open the cassette.

5. Embed, cool and section as per your usual protocol.

6. Dispose of the lid in accordance with laboratory guidelines.

DO NOT REUSE THIS PRODUCT.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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