

CellStor™ Pots

Please
contact us for
more information
cellpath.com

Pre-filled specimen containers

Designed to make the lives of laboratory staff safer and easier, the CellStor™ pot is supplied partially filled with either neutral buffered formalin or formal saline.

Pots are available in 20ml, 60ml, 90ml, and 120ml with respective 13ml, 40ml, 60ml and 90ml fills, come packed into convenient trays, and come complete with relevant material safety data sheet. They are batch numbered and shrink wrapped for extra security.

The shatterproof polypropylene CellStor™ pot contains an innovative sealing mechanism, which eliminates the possibility of leakage during transport, giving laboratories the peace of mind that a sample will reach its destination safely and correctly fixed.



Product Tolerances:

Formaldehyde content to within + or – 0.2%

Specific gravity to within + or – 0.03%

pH to within + or – 0.1% on Neutral Buffered Formalin and + or – 0.75% on Formal Saline

Part No	Description	Qty
BAF-5000-08A	CellStor Pot 20ml - 10% Neutral Buffered Formalin	25
BAF-6000-08A	CellStor Pot 60ml - 10% Neutral Buffered Formalin	25
BAF-9000-08A	CellStor Pot 90ml - 10% Neutral Buffered Formalin	25
BAF-2000-70A	CellStor Pot 120ml - 10% Neutral Buffered Formalin	20
BAF-5100-08A	CellStor Pot 20ml - 10% Formal Saline	25
BAF-6100-08A	CellStor Pot 60ml - 10% Formal Saline	25
BAF-9100-08A	CellStor Pot 90ml - 10% Formal Saline	25
BAF-2100-70A	CellStor Pot 120ml - 10% Formal Saline	20