

## CelliVizor Face Shield

Manufactured from recyclable PET, the multi-use CelliVizor Face Shields bring a number of advantages over masks, such as protecting the eyes, nose and face – which can be an entry point into the body for some viruses. They can also help to reduce the risk of people touching their face and inoculating themselves with any virus they have picked up on their hands.

The smear-free material with antifog coating allows the protective face shield to remain clear throughout its intended use.

Cleaning and sanitising the face shield is easily achieved using soap and water or 70% alcohol wipes. Simply dry using an absorbing cloth to complete the process.

The CelliVizor Face Shield comes ready to use, in a one-piece construction to fit all users.

**To Fit:** Simply place the flexible strap around your head whilst holding the shield to the front of your face. Ensure the shield covers your eyes, nose and mouth. Adjusting the strap to fit is easily done via either of the methods displayed in the 'Instructions for Use' documentation.

Item Code	Description	Quantity
XTB-0100-00A	CelliVizor Face Shield	300



### KEY BENEFITS

- 1 Excellent face and eye protection from splashing or spraying
- 2 Designed to be worn with or without eyewear and face masks
- 3 CE marked; BS EN 166:2002
- 4 Simple to recycle; manufactured from a high-end PET film which is readily recyclable
- 5 Low cost; inexpensive initial outlay and reduced ongoing costs through visors being re-usable
- 6 Easy to clean between use; the complete visor can be cleaned quickly and effectively for re-use in line with current HSE guidelines
- 7 Single piece, optically clear construction; maximises visibility and reduces potential points of contamination
- 8 Space to record individual details; the simple design makes the visor easy to name or label



**Also available**  **Cell Mark**  
- CellMark pens applies a permanent chemically-resistant mark on the face shield.



**CellPath Ltd**  
80 Mochdre Enterprise Park  
Newtown, Powys, SY16 4LE UK  
T: +44 (0)1686 611333  
E: sales@cellpath.com  
[www.cellpath.com](http://www.cellpath.com)