


CellSafe I Biopsy Cassette

The CellSafe I Biopsy Cassette is ideal for the preparation of small biopsies and fragile specimens. Each cassette incorporates an ultra-fine (300 micron), x-ray transparent white mesh that securely accommodates the specimen, whilst permitting efficient processing. This eliminates the need for additional biopsy handling materials such as biopsy pads and bags. The mesh is also acid resistant to permit decalcification.

Manufactured from a specially developed material, the cassette is resistant to all Histology chemicals, retaining its colour during processing. A 45-degree engineered writing surface is optimised for printing on the Leica IPC, Thermo Printmate, and Sakura AutoWrite cassette printers.



Inner chamber rounded corners ensures easy pick-up of biopsies.



Optimal quality printing of text and data matrix barcodes.


The inner chamber has rounded corners throughout, eliminating the possibility of samples 'hiding' in areas hard to get at with forceps.

The rigid structure and design of the cassette significantly increases the adherence of wax blocks when embedding and sectioning.

Dimensions: 28.5mm x 41mm x 7.2mm



LOOSE



HOPPERED



THREADED

The CellSafe I Biopsy Cassette is available either loose, in hopper tubes or threaded.


















KEY BENEFITS

- 1 High specimen security
- 2 Incorporated inner chamber lined with x-ray transparent mesh eliminates the need to open the cassette. This avoids the possibility of losing or breaking the delicate specimen
- 3 Available in a range of colours and packaging options
- 4 Optimised writing surface for printing on the Leica IPC, Thermo PrintMate and Sakura AutoWrite

The CellSafe I Biopsy Cassettes is available in a full range of colours



CellSafe I Biopsy Cassettes

Loose		Hoppered		Threaded		Colour
Item Code	Quantity	Item Code	Quantity	Item Code	Quantity	
EBG-0302-12A	250	EBG-0302-99H	750	EBG-0302-03B	1000	White 
EBG-0304-12A	250	EBG-0304-99H	750	EBG-0304-03B	1000	Green 
EBG-0305-12A	250	EBG-0305-99H	750	EBG-0305-03B	1000	Grey 
EBG-0306-12A	250	EBG-0306-99H	750	EBG-0306-03B	1000	Blue 
EBG-0307-12A	250	EBG-0307-99H	750	EBG-0307-03B	1000	Orange 
EBG-0308-12A	250	EBG-0308-99H	750	EBG-0308-03B	1000	Pink 
EBG-0309-12A	250	EBG-0309-99H	750	EBG-0309-03B	1000	Yellow 
EBG-0310-12A	250	EBG-0310-99H	750	EBG-0310-03B	1000	Red 
EBG-0311-12A	250	EBG-0311-99H	750	EBG-0311-03B	1000	Peach 
EBG-0312-12A	250	EBG-0312-99H	750	EBG-0312-03B	1000	Tan 
EBG-0315-12A	250	EBG-0315-99H	750	EBG-0315-03B	1000	Brown 
EBG-0316-12A	250	EBG-0316-99H	750	EBG-0316-03B	1000	Dark Blue 
EBG-0318-12A	250	EBG-0318-99H	750	EBG-0318-03B	1000	Dark Green 
EBG-0319-12A	250	EBG-0319-99H	750	EBG-0319-03B	1000	Lavender 
EBG-0320-12A	250	EBG-0320-99H	750	EBG-0320-03B	1000	Turquoise 



**CARING FOR OUR
PLANET**

Helping our customers meet their environmental commitments

CellPath are committed to making the world a greener place and keeping our ecological footprint as small as possible.

As such, cassette threads are now manufactured from

specialist cardboard made from pulp that is recyclable and fully biodegradable.

Our hopper tubes are made from rPet (*recycled polyethylene terephthalate*)

which is fully recyclable.



CellPath Ltd

80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: sales@cellpath.com

cellpath.com

