


# Sermo Automated Tape Coverslipper

The Sermo Automated Tape Coverslipper produces coverslipped slides with superior optical quality for reliable long-term storage. The coverslipper is fast, easy to use and efficient, capable of producing approximately 1,200 slides per hour, and is compatible with all commonly used mounting media.

An intuitive user interface features touch-screen controls, allowing up to two slide racks to be loaded simultaneously.



## KEY BENEFITS

- 1 Capable of coverslipping 1,200 slides per hour
- 2 Incorporated drying rack on exit
- 3 2 slide rack capacity
- 4 Adjustable tape length option to adapt to laboratory requirements
- 5 Intuitive software used via touchscreen technology


### Specifications

Item Code:	SAZ-CM1260
Dimensions:	710mm x 750mm x 302mm (d x w x d)

## IPI Pathotape® Coverslipping Tape

**IPI Pathotape® is a high-quality coverslipping film for use in the most common automated tape coverslippers.**

**Possessing high transparency characteristics, Pathotape® provides an excellent impression under the microscope.**



**CellPath**   
INNOVATION IN CELLULAR PATHOLOGY

CellPath Ltd

80 Mochdre Enterprise Park,  
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: sales@cellpath.com

cellpath.com

