

# Tissue Section Debris Brush with Magnetic End

**CellPath's Tissue Section Debris Brush is used to aid microtomy and ensure the clean working area of the microtome when sectioning.**

Through keeping the cassette clamp and knife holder free from debris, section quality is optimised.

The brush incorporates a magnet on the opposite end to the bristle, which can be used as a way of removing the microtome blade from the knife holder. This eliminates the need for workers to use their hands to change the microtome blade, and avoids potential injury within the laboratory.



Item Code	Description	Quantity
JKA-0100-00A	Tissue Section Debris Brush with Magnetic End	1

## KEY BENEFITS

- 1 Ensures a clean working area
- 2 Prevent wax build up – which can cause chattering, and thick and thin sections
- 3 Improved laboratory safety



**CellPath**   
INNOVATION IN CELLULAR PATHOLOGY

CellPath Ltd

80 Mochdre Enterprise Park,  
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: [sales@cellpath.com](mailto:sales@cellpath.com)

[cellpath.com](http://cellpath.com)

