

# IPI Pathotape® Coverslipping Tape


**IPI Pathotape® is a high-quality coverslipping film for use in the most common automated tape coverslippers. Possessing high transparency characteristics, Pathotape® provides an excellent impression under the microscope.**

The innovative surface structure of the tape allows a homogenous distribution of solvent drops and a reduction in the amount of Xylene required for a successful coverslipping. In comparison to other coverslipping films, the amount of Xylene used could be reduced by 10-30% depending on the automated coverslipper type and the usage conditions.

This surface structure also improves long-term adhesion without any effect on the high transparency of the film after the coverslipping procedure.

Rolls are sealed and are 70 meters in length, which allows the coverslipping of approximately 1,400 slides.


IPI Pathotape® is manufactured in Germany and developed and quality controlled according to DIN EN ISO 13485 and under the restrictions of (certified) DIN EN ISO 9001/14001/50001.



## KEY BENEFITS

- 1 Compatible with the most commonly used automated tape coverslippers
- 2 Reduces the amount of Xylene required by 10-30%
- 3 One roll allows the coverslipping of approximately 1,400 slides

Item Code	Description	Quantity
SAZ-0100-00A	IPI Pathotape® Coverslipping Tape – 70m Roll	1



**CellPath**   
INNOVATION IN CELLULAR PATHOLOGY

CellPath Ltd

80 Mochdre Enterprise Park,  
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: sales@cellpath.com

cellpath.com

