

CellWax Plus™

CellWax Plus™ is a highly purified and filtered white translucent paraffin wax that provides great tensile strength and improved elasticity when sectioning.

Our embedding wax is characterised by low resistance to cutting, superb ribboning and no cracking or tearing on floatation.



KEY BENEFITS

- 1 Ultra-filtered and pure to ensure no contamination of specimens or knife damage
- 2 Pelletised to ensure rapid melt-down and convenience of use
- 3 Paraffin wax blended with plastic polymers to produce a high-quality product that improves the elasticity and sectioning properties of each block
- 4 Melting point of 54-57°C – will remain stable when molten for long periods

Item Code	Description	Quantity
GCA-0305-00A	CellWax Plus™ (S) – 5kg	1
GCA-0400-00A	CellWax Plus™ (S) – 4 x 2.5kg	1
GCA-0302-00A	CellWax Plus™ (S) – 20kg Sack	1

CellPath
INNOVATION IN CELLULAR PATHOLOGY

CellPath Ltd

80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK

T: +44 (0)1686 611333

E: sales@cellpath.com

cellpath.com

