





System II Hex Laser Cassettes

Manufactured at our UK headquarters, CellPath's System II Hex Laser Cassette is compatible with the most commonly used laser cassette printers.



Optimal quality printing of text and data matrix barcodes.

The Hex Laser Cassette incorporates a unique and innovative honeycomb design that maximises the open surface area and increases reagent flow through for improved processing. It also results in a more rigid block when sectioning.

The printing face has been optimised for use with a laser cassette printer, resulting in the printing of high contrast text and data matrix codes for near perfect scan rates. Extra thick side walls further increase its stability and rigidity of the cassette, essential for high quality sectioning.

KEY BENEFITS

- 1 Hex design maximises open surface area to increase reagent flow
- 2 Optimised writing surface for printing high-quality text, 1D and 2D barcodes
- 3 Large colour range available



The System II Hex Laser Cassettes are available in a full range of colours



Tissue Processing Cassettes

Easy Loading Threaded Stacks

Our innovative threaded solution manufactured from fully recyclable cardboard allows the easy loading of 40 cassettes into your laser printer via one simple insert and can be used with the majority of available laser printers on the market.



System II Hex Laser Cassettes


Item Code		Quantity	Colour	
EAD-0202-72B	System II Hex Laser Cassettes - White	2000	White	<input type="checkbox"/>
EAD-0204-72B	System II Hex Laser Cassettes - Green	2000	Green	<input type="checkbox"/>
EAD-0205-72B	System II Hex Laser Cassettes - Grey	2000	Grey	<input type="checkbox"/>
EAD-0206-72B	System II Hex Laser Cassettes - Blue	2000	Blue	<input type="checkbox"/>
EAD-0207-72B	System II Hex Laser Cassettes - Orange	2000	Orange	<input type="checkbox"/>
EAD-0208-72B	System II Hex Laser Cassettes - Pink	2000	Pink	<input type="checkbox"/>
EAD-0209-72B	System II Hex Laser Cassettes - Yellow	2000	Yellow	<input type="checkbox"/>
EAD-0211-72B	System II Hex Laser Cassettes - Peach	2000	Peach	<input type="checkbox"/>
EAD-0212-72B	System II Hex Laser Cassettes - Tan	2000	Tan	<input type="checkbox"/>
EAD-0219-72B	System II Hex Laser Cassettes - Lavender	2000	Lavender	<input type="checkbox"/>

System II Hex Cassette Lids

The System II Hex Laser Cassette is used in combination with the System II Hex Cassette Lid, with the choice of an easy to attach disposable plastic lid, or a re-usable stainless steel version depending on the laboratories requirements.

KEY BENEFITS

- 1 Compatible with the System II Hex Cassette
- 2 Two options available; disposable plastic lid, or a reusable stainless steel lid



Item Code	Description	Quantity	Colour	
EAE-0102-03A	System II Hex Cassette Lid - White	1000	White	<input type="checkbox"/>
EAW-0114-08A	System II Hex Cassette Lid - Stainless Steel	25	SST	<input type="checkbox"/>

CellPath 
INNOVATION IN CELLULAR PATHOLOGY



Environment Benefits

CellPath are committed to lowering our ecological footprint, and helping customers achieve the same. Reusable stainless steel lids are therefore offered as an environmentally friendly alternative to plastic lids. Our threaded laser cassette packaging option is also made from recycled cardboard that is fully biodegradable.

CellPath Ltd

80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com

